



Instructions for Transtibial Amputees with Pin Locking Suspension System

Application:

Roll the silicone liner inside out and flatten the end before placing it on your limb. There should be no gapping between the liner and your limb.

Continue to roll the liner onto your limb being careful of the pin (the pin can poke the liner if you are not careful).

The liner should be rolled completely onto the limb and never pulled onto the limb. Tension on the skin may cause skin irritation or blisters. Rolling the top 1/3 of the liner up and down a couple of times will help reduce the tension.

Once the liner is in place continue to add the appropriate amount of prosthetic socks.

- Be sure that the socks are free of any wrinkles.
- Maintaining the correct amount of prosthetic socks throughout the day will be an important part of a successful prosthetic experience.
- If at any time your prosthesis seems too loose or too tight, you will need to modify the amount of socks.
- It is recommended that you carry an assortment of socks with you for changes in volume (size of your limb) throughout the day.

Step into your prosthesis until you are securely locked in

Remove the Prosthesis:

Pushing in and hold the button to release the lock mechanism. If the lock sticks you may have to put weight onto the prosthesis and then disengage the lock.

Remove your limb from the prosthesis.

Remove the liner if applicable.

Remove any prosthetic socks and then roll your liner off of your limb.

If at any time you have questions or concerns, please contact us.